

Protocol Name: **CSF FLOW STUDY**

Body Part: **HEAD**

Coil: **HEAD**

Sequences	Seq Name	Image Plane	# Slices	TR	TE	Slice Thick	FOV	NSA	Matrix	Phase Dir
3 Plane Loc	gre	3 Plane Loc	3	8.6	4	7	250	1	256	A>P
Sag 3D T1	gre	Sagittal	144	19	4.92	1	250	1	256	A>P
Sag Inplane Flow	fl_pc	Sagittal	1	21.3	5.92	5	220	1	256	A>P
Axial Through Plane	fl_pc	Axial	1	26.5	6.9	5	140	2	256	R>L

Details: This Axial Through-Plane series indicates only one slice. The study should be set with 3 Axial slices.  
 Slice One at Foramen Magnum  
 Slice Two spaced evenly between  
 Slice Three at Mid Pons