

Protocol Name: **EOVIST ABDOMEN**

Body Part: **LIVER**

Coil: **BODY ARRAY**

Sequences	Seq Name	Image Plane	# Slices	TR	TE	Slice Thick	FOV	NSA	Matrix	Phase Dir
BH loc	gre	3 Plane	15	7.8	3.69	6	400	1	256	A>P
Axial Inphase/Outphase	fl3d_vibe	Axial	64	4.36	1.33/2.45	3	380	1	320	A>P
Axial T1 FS Pre	fl3d_vibe	Axial	64	3.92	1.39	3	380	1	320	A>P
Eovist Injection										
Axial T1 FS Immediate 1	fl3d_vibe	Axial	64	3.92	1.39	3	380	1	320	A>P
Axial T1 FS Immediate 2	fl3d_vibe	Axial	64	3.92	1.39	3	380	1	320	A>P
Axial T1 FS Immediate 3	fl3d_vibe	Axial	64	3.92	1.39	3	380	1	320	A>P
Axial T1 FS 3 min Post	fl3d_vibe	Axial	64	3.92	1.39	3	380	1	320	A>P
Cor T1 FS 4 min Post	fl3d_vibe	Coronal	128	3.16	1.12	1.8	400	1	288	R>L
Cor T2 BH	haste	Coronal	30	1400	93	5	380	1	256	R>L
Axial T2 BH	haste	Axial	30	1400	93	5	380	1	256	A>P
Axial T2 FS BH	haste	Axial	35	2000	93	5	380	1	320	A>P
Axia Diffusion	ep2d_diff	Axial	25	5900	83	5	380	5	192	A>P
Axial T1 FS 20min Post	fl3d_vibe	Axial	64	3.92	1.39	3	380	1	320	A>P
Cor T1 FS 20min Post	fl3d_vibe	Coronal	128	3.16	1.12	1.8	400	1	288	R>L

Details: Diffusion - B-values are 50,400,1000