

Protocol Name: **FOREFOOT WITH AND WITHOUT CONTRAST**

Body Part: **FOOT**

Coil: **MEDIUM FLEX COIL**

| Sequences | Seq Name | Image Plane | TR | TE | Slice Thickness | FOV | NSA | Matrix | Phase Dir |
|-------------------|----------|-------------|----|----|-----------------|-----|-----|--------|-----------|
| Axial Localizer | GRE | Axial | | | | | | | |
| 3 Plane Localizer | GRE | 3 Plane | | | | | | | |
| Axial T1 | FSE | Axial | | | | | | | |
| Axial T1 FS Pre | FSE | Axial | | | | | | | |
| Axial T2 FS | FSE | Axial | | | | | | | |
| Sag T1 | FSE | Sagittal | | | | | | | |
| Sag T2 FS | FSE | Sagittal | | | | | | | |
| Cor T2 FS | FSE | Coronal | | | | | | | |
| Cor T1 | FSE | Coronal | | | | | | | |
| Gadolinium | | | | | | | | | |
| Axial T1 FS Post | | Axial | | | | | | | |
| Sag T1 FS Post | | Sagittal | | | | | | | |
| Cor T1 FS Post | | Coronal | | | | | | | |

Details:

1. Axials are oriented axially to forefoot
2. Coronals are oriented Plantar to Dorsal