

Protocol Name: **HUMERUS ROUTINE**

Body Part: **ARM**

Coil: **BODY FLEX COIL**

Swquences	Seq Name	Image Plane	# Slices	TR	TE	Slice Thick	FOV	NSA	Matrix	Phase Dir
Localizer	gre	2 Plane loc	3	7.8	3.69	6	450	1	256	R>L
3 Plane Loc	gre	3 Plane loc	9	7.8	3.69	6	360	1	256	A>P
Cor PD	tse	Coronal	25	3000	28	3.5	360	1	320	R>L
Cor T2 FS	tse	Coronal	25	4000	47	3.5	360	1	320	R>L
Sag T1	tse	Sagittal	30	763	17	3	360	2	320	A>P
Sag T2 FS	tse	Sagittal	30	3500	59	3	360	2	320	A>P
Axial T2 FS Upper	tse	Axial	30	5270	100	4.5	160	1	384	A>P
Axial T2 FS Lower	tse	Axial	30	5270	100	4.5	160	1	384	A>P
Axial T1 upper	tse	Axial	30	550	21	4.5	160	1	320	A>P
Axial T1 Lower	tse	Axial	30	550	21	4.5	160	1	320	A>P