GE VCT 64 SLICE PROTOCOLS

Protocol Name:]	INVASIVE PROCEDURE ADULT (Select routine protocol for body part ordered)												
NDICATIONS:			Injection of contrast into joint for MRI or steroid injection into joint				LABEL:				Sa	Same as in protocol selected.			
IV CONTRAST:			Usually none				RECON 1:				5mm thick / 5 mm interval. Storecon 400ww/30wl or ask radiologist				
ORAL CONTRAST:			Usually none.				RECON 2 AND ADDITIONAL POST PROCESSING:				No				
PATIENT POSITIONING AND SCAN AREA:			Position patient as instructed by radiologist. Select routine study for the body part. Take your scout and images through the area of interest with grid marker. Select biopsy Rx once centered on the grid marker. Mark the best spot for the needed with a sharpie once Rad tells you where to mark.				THINGS TO REMEMBER:				Need procedure consent form signed by patient after radiologist explains procedure and site verification form with "time out" time recorded for an needle stick procedure.				
GENERAL EXAM			Usually 5 x 5 mm axial images-				SEND EXAM:				Send scout, axial images, and procedure images to pacs.				
GUIDELINES:			ask radiologist. DFOV large enough to see outside of skin.								pro	ocedure	images to	o pacs.	
RIGGER		F A 3.7	Usually none				ARCHIVE:								
CONTRAST DELAY: FECHNIQUE GUIDELINES:			Varies by protocol selected				BILL AS:				Varies. Charge for procedure, sterile tray, needles used, & sedation if needed.				
				SCOU	T PRC	TOCO	L SET	IN S	CANNE	ER					
Scan Type		Start Loc		End Loc		kV		mA			out		Voic Light Time	s	
Varies															
				XIAL IN											
Scan Type	Start Loc	End Loc	Thick Speed	Interval (mm)	Gantry Tilt	SFOV	kv	ma	DFOV	Show Recon1		Show Rocon2	Show Rocon3	Direct Vis	
Varies															