

Protocol Name: **PELVIS ROUTINE**

Body Part: **PELVIS**

Coil: **BODY ARRAY**

Sequences	Seq Name	Image Plane	# Slices	TR	TE	Slice Thick	FOV	NSA	Matrix	Phase Dir
Localizer	gre	3 Plane	7	7.8	3.69	6	400	1	256	A>P
Sag T2	tse	Sagittal	48	5960	94	5	300	2	384	H>F
Sag T1	tse	Sagittal	48	650	11	5	300	1	384	H>F
Cor T2 FS	tse	Coronal	40	6680	55	4	340	2	384	R>L
Cor T1	tse	Coronal	40	726	11	4	340	2	384	R>L
Cor GRE	gre	Coronal	40	970	19.9	4	340	1	256	R>L
Axial T1	tse	Axial	48	775	11	4	340	2	384	A>P
Axial T2 FS	tse	Axial	48	4140	55	4	340	2	384	R>L

