

3T MRI PROTOCOLS

Protocol Name: **PROSTATE W/WO CONTRAST**

Body Part: **PELVIS**

Coil: **BODY ARRAY**

Sequences	Seq Name	Image Plane	# Slices	TR	TE	Slice Thick	FOV	NSA	Matrix	Phase Dir
Localizer breath-held	gre	3 Plane	7	7.8	3.69	6	400	1	256	A>P
COR T2 BH	haste	Coronal	44	1200	122	5	440	1	320	R>L
Axial T2 BH	haste	Axial	64	1200	122	5	380	1	320	A>P
Localizer breath-held	gre	3 Plane	7	7.8	3.69	6	400	1	256	A>P
Cor T2 SPACE FS	tse	Coronal	192	1600	108	1	380	1	384	R>L
Sag T2	tse	Sagittal	35	4840	108	3.5	160	2	320	R>L
Axial T2	tse	Axial	40	5400	108	3	160	1	320	R>L
COR T2	tse	Coronal	35	4850	108	3.5	160	1	320	R>L
Axial T1 PRE	tse	Axial	40	505	12	3	160	1	192	R>L
Axial T2 FS	tse	Axial	48	4140	55	4	340	2	384	R>L
R>L Axial DWI	DWI	Axial	40	7600	90	3	200	3	102	R>L
Axial T1 VIBE test	VIBE	Axial	24	5.07	1.79	3.5	160	1	192	R>L
Gadolinium										
Axial T1 VIBE DYN 35passes	VIBE	Axial	24	5.07	1.79	3.5	160	1	192	R>L
Axial T1 FS Post	tse	Axial	40	635	12	3	160	1	192	R>L