

3T MRI PROTOCOLS

Protocol Name: **RENALS NON-CON NATIVE**

Body Part: **ABDOMEN**

Coil: **BODY ARRAY**

Sequences	Seq Name	Image Plane	# Slices	TR	TE	Slice Thick	FOV	NSA	Matrix	Phase Dir
Localizer BH	gre	3 Plane	15	7.8	3.69	6	400	1	256	A>P
Cor T2 BH	haste	Coronal	30	1400	93	5	380	1	256	R>L
Axial T2 BH	haste	Axial	30	1400	93	5	380	1	256	A>P
3 Plane Loc	gre	3 Plane	33	6.6	3.69	6	460	1	256	A>P
3 Plane loc for Navigator	gre	3 Plane	21	7.8	3.69	6	400	1	256	A>P
NATIVE_TrueFISP_nav_ECG	CV_nav	Axial	88	750.65	1.7	0.9	340	1	304	A>P
NATIVE_TrueFISP_nav_NO ECG	CV_nav	Axial	72	1277	1.55	1	320	1	304	A>P

Details: