

3T MRI PROTOCOLS

Protocol Name: **SHOULDER ARTHROGRAM**

Body Part: **SHOULDER**

Coil: **8 CHANNEL**

Sequences	Seq Name	Image Plane	# Slice	TR	TE	Slice Thick	FOV	NSA	Matrix	Phase Dir
2 Plane Localizer	gre	2 Plane	6	7.8	3.69	6	280	1	256	R>L
Axial T2 FS	tse	Axial	25	4000	53	3	140	1	256	A>P
Axial T1 FS	gre	Axial	25	761	11	3	140	1	256	A>P
Cor Obl T1	tse	Coronal Obl	20	600	22	3	140	1	384	F>H
Cor Obl T1 FS	tse	Coronal Obl	20	608	19	3	140	1	256	F>H
Cor Obl T2 FS	tse	Coronal Obl	25	4000	48	3	140	1	256	F>H
Sag Obl T1 FS	tse	Sagittal Obl	20	657	19	3	140	1	256	F>H
ABER										
2 Plane Localizer	gre	2 Plane	6	7.8	3.69	6	280	1	256	R>L
ABER Sag Obl T1 FS	tse	Sagittal Obl	20	657	10	3	140	1	256	F>H