

Protocol Name: **WRIST ROUTINE**

Body Part: **WRIST**

Coil: **8CH WRIST COIL**

Sequences	Seq Name	Image Plane	# Slices	TR	TE	Slice Thick	FOV	NSA	Matrix	Phase Dir
Axial Localizer	gre	Axial	5	7.8	3.69	6	400	1	256	A>P
3 Plane Localizer	gre	3 Plane	9	8.3	3.69	6	200	1	256	A>P
Axial T2 FS	tse	Axial	25	4100	52	2	70	1	320	A>P
Axial T1	tse	Axial	25	700	15	2	70	1	384	A>P
Cor DESS	dess	Coronal	104	16	5.2	0.4	100	1	256	R>L
Cor T1	tse	Coronal	16	520	30	3	70	1	384	R>L
Cor PD FS	tse	Coronal	16	3000	39	3	70	1	320	R>L
Sag T2 FS	tse	Sagittal	20	4100	77	3	90	1	320	A>P

Details: Patient in "Superman" Position with wrist in coil, arm of interest extended over head